

Adult Crisis Auscultation Manikin
Item# U56930FLALD

This manikin has all the same features as the complete Crisis Manikin and also includes auscultation sites with heart and lung sounds. For auscultation training, the instructor uses the wireless remote to select from a menu heart and lung conditions. The manikin appears as a real patient without any visible auscultation sites. The student must palpate to identify correct auscultation locations and will hear different heart and lung sounds. A diagnosis of the condition selected by the instructor can be made by comparing the variations in sounds occurring at different sites, just like a real patient. Lung sounds can be detected at five anterior and two mid-axillary locations and students can practice auscultation at six anterior heart sites. The remote control does not have to be pointed directly at the manikin or to operate.

Features:

The manikin includes Complete Crisis Manikin with auscultation sites, one remote control with LCD display, one scope with single- and dual-user headpieces, and hard storage case. Requires two AA and two AAA batteries (included). Ideal for ACLS, paramedic, EMT, and nursing training at every level. Five-year warranty on the manikin.

Lung Conditions:

- Bronchial
- Egophony
- Bronchovesicular
- Fine Crackle
- Cavernous
- Friction Rub
- Coarse Crackle
- Infant
- Mono Wheeze
- Normal
- Normal Vesicular
- Pectoriloquy
- Pulmonary Edema
- Rhonchi
- Stridor
- Wheeze



Heart Conditions:

- Normal
- Aortic Regurgitation
- Mitral Stenosis
- Pulmonary Stenosis
- Holosystolic
- Mid Systolic
- S3 and S4 Gallups
- Systolic Click
- Atrial Septal Defect
- PDA & VSD

