

**Bariatric Patient Lift 1000 LB Capacity, Crank  
Item# 0001PTUF**

Patient Lift is able to lift bariatric patients and provides them mobility without the use of an expensive hoist system.

**Features:**

- The frame is constructed out of a strong carbon steel with its high stress areas reinforced.
- It uses a durable heavy duty hydraulic cylinder that is operated by a crank.
- Features a lifter arm with adjustable length to adapt to your transfer needs.
- Lockable rear casters ensure operation stability and keep the 0001P sturdily in place.
- Foot pedal allows caregiver to manipulate the base legs without needing to bend over and adjust the base manually.



**Specifications:**

- Pivot height x overall length: 58in x 63in
- Base size: 53in x 36in
- Minimum underbed clearance: 4.5in
- Minimum ceiling height: 86in



Toll Free 1.877.706.4480 | Local/ International 713.706.4480 | Fax 713.706.4161  
[www.MSEC.co](http://www.MSEC.co)