



toll free: 1-877-706-4480
international/local: 713-706-4480
fax: 713-706-4161
www.MS3C.co

Battery Powered Stair Climb Assist Chair
Item# MS3C-300TSB

MS3C-300TSB Battery Powered Stair Evacuation Chair is ideal for emergency evacuation during fires or power outages or for use at home for assisted mobility up and down the stairs. Electric powered unit significantly minimizes caregiver injury by replacing manual lifting with electric operation. With high durability rubberized track, individuals can be easily transported up or down on virtually any staircase surface. Ideal for medical facilities, multi-story commercial and residential buildings, fire and EMS vehicles.

Features:

- Designed for fast and safe stair evacuation in both residential and commercial settings
- Light-weight, aluminum construction powder-coated frame
- Easy to read color coded navigational controls with night light
- Ergonomic lifting handles with molded hand grips to reduce stress and provide extra control for caregiver
- Locking guide handles extend for secure lifting if necessary
- Telescoping lift handles easily adjust to caregivers specific needs
- Dual, extra traction treads withstand high volume use, high durability tracks for non-slip traction on any surface with adjustable track
- 4 in swiveling front castors increase range of motion and sharper turn radius for use in confined spaces
- Dual locking rear wheels with sealed bearings stabilize patient loading and unloading
- Easy-to-fold, compact for convenient storage for small spaces.
- Retractable solid foot rest ensures safe transport and added comfort
- 4-point patient harness included
- Included head support and headrest is easily removed for storage
- Easy to clean, grease-free maintenance



Included Accessories:

- Color Coded buttons for easy navigation
- Night light
- Neoprene seat top on a solid platform for extra comfort
- Neoprene back cushion



toll free: 1-877-706-4480
international/local: 713-706-4480
fax: 713-706-4161
www.MS3C.co

- Four point patient safety restraints
- Head rest strap
- Head rest cushion
- Battery charger and power cord
- Heavy Duty Lithium-Ion battery

Specifications:

- Minimum recommended stair case landing 42 x 42in
- Folded size (inches): 45L x 19.5W x 10.6H in
- Net weight: 66 lbs
- Static weight capacity: 500 lbs
- Operational weight capacity: up to 400 lbs
- Battery operating time: 1.5-2 hours (Equivalent to 200+ flights of stairs)
- Motor: 200W 24V 13.9A
- Output: 20-29.4V
- Replacement batteries sold separately
- FDA, CE and ISO 9001 Certificate



Limited Manufacturer's Warranty: 3 year on frame, 1 year on motor and controller

Accessories Available:

- Dust Cover
- Wall mount brackets
- EZ Store

