

Chiropractic 2-Section Hi-Lo Bariatric Table

Features:

- Extra-wide 40in base and cushion
- Gas-assisted struts control the recessed headrest
- Table height electrically adjusts from 18in-40in for ease of access
- Foot-pedal control on all four sides standard
- Industry-leading, linear-actuator silently generates 7,000 Newton force
- 110V or 220V options with the flip of a switch
- The 3in thick non-slip cushion is gently contoured for both patient and practitioner comfort
- Integrated, retractable casters included
- Available in black, gray or blue (other colors available for an additional fee)




Specifications:


- Height: 18 in to 40 in
- Width: 40 in
- Length: 80.7 in total Head: 25.6 in, Leg 55.1 in

Section Configuration: (see diagram below)

- Head: 25.6 in
- Leg 55.1 in



Designed and Engineered in
the United States



MEDICAL EQUIPMENT & CONSTRUCTION SPECIALISTS

21830 Kingsland Blvd-Suite 101A Katy TX 77450

P 713.706.4480 F 713.706.4161

www.msec.co