

# Chiropractic Hi-Lo Table with Stool


## Features:

- Gas-assisted struts control
- Four cushion-sections
- Low-friction gliding lumbar section
- Turret-mounted traction pedestal
- Table height electrically adjusts from 18in-38in for ease of access
- Foot-pedal control on all four sides
- 110V or 220V are options with the flip of a switch
- Antimicrobial surface provides a sanitary surface that is easy to clean
- 3 in non-slip cushion is gently contoured for both patient and practitioner comfort
- Integrated retractable casters included
- Available in black, gray or blue




## Specifications:

- Height: 18 in to 38 in
- Width: 27.2 in
- Length: 81.4 in
- 4 sections: 15.6 in, 11.9 in, 11.9 in, 28.3 in
- Step stool: 19.7 in x 11.8 in



## Section Configuration: *(see diagram below)*



21830 Kingsland Blvd-Suite 101A Katy TX 77450  
P 713.706.4480 F 713.706.4161  
www.msec.co

