

# Aluminum Alloy Ambulance Stretcher

Model: MS3C-200C



This EMS/ Ambulance gurney is designed for both commercial and private use. Durable enough to withstand the rigors of commercial patient transportation while easy enough to operate by a single individual.

## Standard Features:

- Made of aluminum alloy and stainless steel materials
- Stretcher comes standard with heavy duty, oversized locking casters
- Mechanical folding legs are controlled by the handles on both sides
- Siderails collapse allowing for easy patient loading and unloading
- Only one operator needed to load the patient into the emergency vehicle
- 2 inch waterproof foamed cushion is adjustable
- Light-weighted, easy to sterilize and clean
- Pneumatic backrest and footrest allows for easy adjustments
- Stretcher release handles are conveniently placed for easy operation
- Floor-mounted vehicle locking device locks stretcher while in motion
- Reinforced bars protect the stretchers legs when loading and unloading the patient
- (2) Patient safety belts included



## Specifications:

- Max vehicle loading height: 27in
- Upright dimensions (w/ pad): 77.5in x 24in x 37in
- Collapsed dimensions (w/ pad): 78in x 24in x 13in
- Packing size: 82in x 24in x 19in
- Load capacity: 400 lbs.

**Manufacturer warranty:** 1 year

**Note:** Stretcher's O2 holder ships inverted. To use, remove screws and re-install with holder facing out. Shown below, shipped and installed.

Shipped Inverted



Flipped and Installed

