



toll free: 1-877-706-4480
international/local: 713-706-4480
fax: 713-706-4161
www.MS3C.co

Lite Emergency Crash Cart with Accessory Package and IV Pole

Item# MS3C-500CR

The emergency cart makes organization easy. Standard features and accessories provide dependable performance and promote improved workflow in emergency departments. The breakaway locking handle offers security while still providing easy access necessary to every emergency response. Stable maneuverability is built-in to every cart and comes standard with a full line of accessories enhance user efficiency and support code team performance.

Features:

- Made of ABS plastic, easy to clean and corrosion-proof
- All-in-one table top with raised sides to eliminate items rolling off
- Built-in handrail for easy maneuvering
- ABS table top with transparent dust-proof pad
- Sliding side shelf for writing
- 5 ABS drawers slide smooth and silent with centralized lock and plug-in labels to avoid dust
- (3)Flexible material storage bins that can be used for gloves, etc.
- (4) Anti-Static castors, two locking, designed to eliminate debris from wrapping the caster
- Drawer dividers for arbitrary partitioning
- (2)Trash bins
- Telescoping IV pole
- Defibrillator Holder
- Power Cord
- C.P.R board
- Oxygen tank holder

Drawer Configuration:

- (2) Small drawers: 20in (W) x 14in (D) x 2.25in (H)
- (2) Medium drawers: 20in (W) x 14in (D) x 4.5in (H)
- Large drawer: 20in (W) x 14in (D) x 8.25in (H)

Dimensions:

- 19 in (D) x30 (W) in x 38in (H)

