



toll free: 1-877-706-4480
international/local: 713-706-4480
fax: 713-706-4161
www.MS3C.co

Manual Hydraulic Surgical Table with Surgical Accessories Item# MS3C-800HST

This state of the art universal operating table combines the latest technology with the most affordable pricing in the market.

Features:

- Manually operated functions for optimal performance
- Stainless steel columns that provide stability
- Modern sleek look for any surgical suite
- The versatility allows the surgeon to perform many types of surgical procedures such as gynecological, obstetrics, ophthalmology, orthopedic, urology, abdominal and chest surgery
- Hydraulic pump system elevates the table with a foot pedal making it extremely convenient as well as more reliable than electrically operated tables
- All Trendelenburg functions and lateral movements are controlled by a mechanical operating system which can be manually manipulated with handles conveniently located on the table's midsection
- 2in thick high-density mattress
- CT scan Compatible



Standard Accessories:(included):

- 1 Head rest
- 2 Shoulder supports
- 2 Arm boards
- 2 Lateral supports
- 2 Stirrups/Knee Crutches

Specifications:

- Tabletop length and width: 82in x 19in
- Elevation: 31-41in
- Trend/ Reverse trend: +25 degree
- Lateral tilt: +20 degree
- Leg section adjustment: +90 degree
- Head section adjustment: 30 degree / 90 degree
- Back Plate: 90 degree / 10 degree

Manufacturer warranty: 1 year