



toll free: 1-877-706-4480  
international/local: 713-706-4480  
fax: 713-706-4161  
[www.MS3C.co](http://www.MS3C.co)

## **Pediatric Tilt Wheelchair**

### **Item# MS3C-PTWC**

#### **Features:**

- Newly designed reinforced frame
- Newly designed height adjustable & removable "T" style armrests
- Folding frame tilt-in-space wheelchair with 45 degrees of tilt
- Removable anti tippers included
- Hand Brake included
- Transit brackets included for securing unoccupied wheelchair in vehicle
- Included accessories: Head rest, Foot rests, Lateral Trunk support, Removable anti tippers included, Thigh and Trunk Support, Patient Safety belt
- Specifications:
- Weight capacity: 180 lbs
- Seat width: 15 in
- Seat Depth: 13 in
- Locking Rear Wheels: 15 in
- Front wheels: 6 in
- Weight: 45 lbs with accessories
- Weight : 28 lbs without accessories



#### **Dimensions:**

- 37 in H x 19 in W x 38 in D

#### **Folded Dimensions:**

- 37 in H x 11 in W x 34 in D