

## Tall Slim Aluminum Alloy Crash Cart with Accessories

### Item# MS3C-610CRM

The emergency cart makes organization easy. Standard features and accessories provide dependable performance and promote improved workflow in emergency departments. Stable maneuverability is built-in to every cart and comes standard with a full line of accessories enhance user efficiency and support code team performance.

#### Features:

- Four anti-scratch aluminum alloy columns with ABS body
- ABS plastic top board with stainless steel side rails around for protection
- IV Pole with 4 I.V. transfusion hook at the top
- With disposable sharp container, trash bin, extending sliding writing shelf, utility container and wire basket
- Four 4 in castors and two with locking brakes
- 5 Drawers with dividers
- O2 Holder and power supply included
- Color: Orange

#### Cart Dimensions:

- 41 in H x 19 in D x 26 in W

#### Drawer Configuration:

- 5 drawers: 3(small) 1(medium) 1(large)

#### Drawer Specifications:

- 4in H x 14in D x 19in W (small)
- 5.5in H x 14in D x 19 holder19in W (medium)
- 7in H x 14in D x 19 holder19in W (large)

**Manufacturer warranty: 1 year**

